


CONSUMABLES

Reagent Bags

The Pool Lab ASP reagent bags contain 150mL of reagent each, which will last for approximately 560 tests, or approximately 6 months at 3 tests per day. If the bags are empty, then the ASP cannot perform water tests. In this case there will be no acid dosing, and the chlorinator will continue to run based on previous demand.

NOTE: The Free Chlorine Reagent #1 has a limited shelf life and will slowly change colour over time. Do not be concerned if this reagent has a magenta colour when you install it. The reagent is considered expired only when it begins to turn brown and is no longer transparent.

RED = pH REAGENT **WHITE = CHLORINE REAGENT #1** **BLUE = CHLORINE REAGENT #2**



NOTE: After installing new reagent bags, prime the reagent lines by either:

1. Cycle power OFF and back ON after changing reagent bags (clears test history)
2. Select PRIME REAGENT function from the menu system.
 - (a) Go to the MAIN MENU (press and hold VIEW for 3 seconds)
 - (b) Navigate to ASP TOOLS, press VIEW
 - (c) Navigate to REAGENTS, press VIEW
 - (d) Select PRIME REAGENTS, press VIEW
 - (e) Press VIEW once more to initiate the prime sequence.